

Human Saliva, Toxicology Neg.

Catalog Number: OHI060-DF

Lot Number: <Lot Specific>

Storage: -20 °C
Manufacture Date: <Lot Specific>

Source: Human Saliva

Product Performance Specifications

The following analytes screened by LC-MS/MS, results are None Detected:

d-methamphetamine	Nordiazepam	Benzoylcegonine	Temazepam	Hydrocodone
d-amphetamine	Delta-9-THC	Morphine	Oxycodone	Hydromorphone
l-methamphetamine	Delta-9-THC-COOH	Oxazepam	Oxymorphone	Cocaethylene
l-amphetamine	11-Hydroxy Delta-9 THC	6-MAM - Free	Chlordiazepoxide	Cocaine
Codeine	Dihydrocodeine / Hydrocodol	Alprazolam	Clonazepam	Dextromethorphan
Diazepam	Dihydrocodeine - Free	EDDP	Lorazepam	Hydrocodone - Free
MDMA	Midazolam	Oxycodone - Free	Phencyclidine	Hydromorphone - Free
MDA	Morphine - Free	Oxymorphone - Free	Methadone	N/A

The following analytes screened by GCMS, results have None Detected:

Amobarbital	Methaqualone	Acetaminophen	MDMA	Nicotine
Methadone	MDEA	Pentobarbital	Propoxyphene	Butabarbital
MDA	Phenobarbital	EDDP	Fentanyl	Fluoxetine
Butalbital	Phencyclidine	Secobarbital	Cotinine	Fluvoxamine
10-Hydroxycarbazepine	Bupropion Metabolite	Desipramine	Doxepin	Gabapentin Breakdown Product
Acetyl Fentanyl	Bupirone	Desmethylcitalopram	Doxylamine	Guaifenesin
Alfentanil	Butorphanol	Desmethylclomipramine	Duloxetine	Haloperidol
alpha-PVP	BZP	Desmethyldoxepin	Duloxetine Breakdown Product	Hexobarbital
Amitriptyline	Carbamazepine	Desmethylmirtazapine	Ephedrine / Pseudoephedrine	Hydroxybupropion
Amlodipine	Carisoprodol	Desmethylsertraline	Eszopiclone / Zopiclone Breakdown	Hydroxychloroquine
Amoxapine	Chlorpromazine	Desmethyltrimipramine	Ethinamate	Hydroxyzine
Amphetamine	Citalopram / Escitalopram	Dextro / Levo Methorphan	Ethosuximide	Ibuprofen



Golden West Diagnostics, LLC · 42191 Sarah Way, Temecula, CA 92590 · +1-951-699-0957

Human Saliva, Toxicology Neg.

Catalog Number: OHI060-DF

Lot Number: <Lot Specific>

Storage: -20 °C

Source: Human Saliva

Manufacture Date: <Lot Specific>

The following analytes screened by GCMS, results have None Detected:

Atomoxetine	Clomipramine	Dextrorphan / Levorphanol	Ethotoin	Ketamine
Atropine	Clonidine	Dicyclomine	Ethylmorphine	Lacosamide
Benzotropine	Clozapine	Diltiazem	Felbamate	Lamotrigine
Brompheniramine	Cyclobenzaprine	Diphenhydramine	Flecainide	Laudanosine
Bupropion	Descarboethoxylorata dine	Donepezil	Fluconazole	Thioridazine
Levamisole	Methylphenidate	Norketamine	Procainamide	Thiothixene
Levetiracetam	Metoclopramide	Normeperidine	Promazine	Ticlopidine
Lidocaine	Metoprolol	Norpseudoephedrine / Phenylpropanolamine	Promethazine	Topiramate Breakdown Product
Loratadine	Metronidazole	Nortriptyline	Propofol	Tramadol
Loxapine	Mexiletine	O-Desmethylvenlafaxine	Protriptyline	Tramadol Metabolite
mCPP	Milnacipran	Olanzapine	Pyrilamine	Trazodone
Meclizine	Mirtazapine	Ondansetron	Quetiapine	Trihexyphenidyl
Memantine	N-Acetylprocainamide	Papaverine	Quinidine	Trimethoprim
Meperidine	Naltrexone - Free	Paroxetine	Quinine	Trimipramine
Mepivacaine	Naproxen	Pentazocine	Rufinamide	Venlafaxine
Meprobamate	Naproxen Breakdown Product	Pheniramine	Sertraline	Venlafaxine Metabolite
Mescaline	Norclozapine	Phensuximide	Strychnine	Verapamil
Mesoridazine	Norcodeine	Phentermine	Sufentanil	Xylazine
Metaxalone	Norcyclobenzaprine	Phenylethylmalonamide (PEMA)	Tetrahydrozoline	Yohimbine
Methamphetamine	Norfentanyl	Phenyltoloxamine	TFMPP	Zaleplon
Methocarbamol	Norfluoxetine	Phenytoin	Theophylline	Zolpidem
Methoxetamine	Norhydroxyzine	Primidone	Imipramine	Zonisamide



Golden West Diagnostics, LLC · 42191 Sarah Way, Temecula, CA 92590 · +1-951-699-0957

Human Saliva, Toxicology Neg.

Catalog Number: OHI060-DF

Lot Number: <Lot Specific>

Storage: -20 °C

Source: Human Saliva

Manufacture Date: <Lot Specific>

The following analytes screened by EIA, results are None Detected:

Benzodiazepines	Cannabinoids	Opiates	Oxymorphone
Buprenorphine / Metabolite	Cocaine / Metabolites	Oxycodone	N/A

Test	Specifications	Methodology
Salicylates	None Detected	Colorimetry
Ethanol	None Detected	GC

Testing

Donors have been tested and found negative or non-reactive for HIV 1/2, HBsAg, HCV, HIV-1 (NAT), HCV (NAT) and RPR by FDA Approved Methods.

Caution

In-vitro use only. Products are intended for research or manufacturing of non-injectable products only. This material was tested as stated above in the testing section. Nevertheless, products of human origin should be considered potentially infectious and handled according to universal precautions. The safe and responsible handling of the product listed above is the responsibility of the customer. We accept no liability for injury or illness resulting from the use and/or handling of the product. The customer assumes responsibility for the proper disposal of all materials in accordance with Federal and State regulations.

Warranty

Golden West Diagnostics, LLC makes no claim or warranty that this product will perform or is suitable for use in a customer's application. It is the customer's responsibility to evaluate and qualify the product's suitability for the customer's intended use or application. Furthermore, this product is a biological material and lots may vary so it is the customer's responsibility to determine if an evaluation of each lot is warranted.



Golden West Diagnostics, LLC · 42191 Sarah Way, Temecula, CA 92590 · +1-951-699-0957